











Medical Bedside Table Instruction



◆ PARTS

| | | | |
|---|---|---|--|
| ①  | ②  | ③  | ④  |
| ⑤  | ⑥  | ⑦  | ⑧  |
| ⑨  | ⑩  | | |

◆ ASSEMBLY

STEP 1

Line up the central hole on the horizontal bar with the bottom hole on the L-shaped black bar.
Put the bolt into the hole, and screw in the nut and tighten it with your hand or wrench.
Install the other horizontal bar in the same way.



STEP 2

Put the threaded head of the wheel into the side hole of the horizontal bar and screw it in with your hand or wrench. And then install the rest 3 wheels in the same way.



STEP 3

Insert the L-shaped silver bar into the L-shaped black bar.

Put the threaded head of the knob into the hole on the black bar and screw it in to tighten the silver bar in place.



STEP 4

Put the slim silver tube into the horizontal tube and press it with one of your hands so that it can stick out of the other end.

And then insert the stop collar into the slim tube and slightly tighten the knob.



STEP 5

Install the tabletop to the horizontal tube. Adjust the angle of tabletop to your desired position and fully tighten the side knob. Then the table is ready for use.



◆ **ADJUSTMENT INSTRUCTION**



1. To adjust the angle of the tabletop, slightly loosen the knob and adjust the tabletop to your desired position, and then tighten the knob.
2. To adjust the height of the table, hold the silver bar and slightly unscrew the knob. Lift or lower the silver bar to your desired position and then tighten the knob.