

2026.5.25

140*210mm

Dear users:

Thank you for choosing our Smart Body Analysis Scale.

This product is designed by your body relationship, to monitor the fat, water, muscle, bone, BMI, visceral fat and so on. The product use 2.4G HZ wireless transmission protocol, can connect Bluetooth from mobile or ipad to transmit data.

Product Features

- 1.Measure Weight, Body Fat, Hydration, Visceral Fat, Muscle Mass, Bone Mass, Calories , BMI and so on.
- 2.Records data to App allowing for you to check progress and maintain goals.

More Details

- Weighting scope/dividing value:5kg-180kg/0.05kg(5-100kg)/0.1kg(100-180kg)
- Battery:2XAAA battery
- Three units: kg/st/lb
- Overload display"ERR", low battery display"Lo"

How to Ensure Accurate Readings

- Always measure barefoot. Do not step on the scale with wet feet, wearing socks or shoes to prevent slipping, damage and inaccurate readings.
- Recommended to be placed on hard and flat surface only.
- Always weigh yourself on the same scale, same place and same surface if possible.
- Recommended measurement time should be after waking up but before the breakfast, and next weight should be two hours after your meal but before your next meal, or two hours after dinner but before bed.
- Erroneous measurements can occur if feet are dirty, dry or extremely callous.
- Body fat percentage estimates will vary with the amount of water in the body, and can be affected by dehydration or over-hydration due to such factors as alcohol consumption, menstruation, illness, intense exercise.
- Stand upright and not to the edge of the scale to ensure good contact with the film or electrodes (the film and electrodes both are conductive media depending scale models).
- Horizontal line should be running below arch of each foot.
- Full feature of scale only available with corresponding app.

Maintenance

- 1.Keep scale surface dry to avoid accident by slipping..
- 2.Please put the scale in a dry place to protect the internal electronic components.
- 3.Please use micro cloth to clean the scale, don't immersed in water.
- 4.Please don't remove the scale except replace the battery.
- 5.Please don't use heavy percussion or falls from a height.

HOW TO USE

1. Installing the battery

1.1. Remove the battery cover at the back of Smart Body Analysis Scale, and insert 2 x1.5V AAA batteries (Not included) into the unit. Make sure the polarity is correct when done.

1.2. To set "units", press the button located on the back of the scale. Choose the right unit you need, kg , st or st:lb. Then, put back the battery cover.

2. Installing app

2.1 Switch on Bluetooth

2.2 Phone or Tablet connects to Net (WiFi,4G or 3G)

2.3 Use your device to scan QR code or search "OKOK International" from

App Store or Google Play to install App.

Click "Turn All Categories On" during installation and allow "okok international" to access all data.



3. Pair the devices and register , complete user information

1. Open "OKOK International"

2. Enter your email Click "Join Us", then input Verification Code and set your password, complete your registration information step by step

Completing a personal profile is essential because the data you enter will be used as reference for APP to analyze and assess health situation.

- If you have already registered before, you can login to your account
- If you are registered and forgot the Password, please click "problems?" and select "Retrieve password". Type your Email, a verification code will be sent to your email. Type the code obtained from your e-mail and reset password step by step.

4. Connect and use weighing scale

- After registration you can start to use your scale to weigh, Click "Add device" and select "weigher", then stand on the scale, after some seconds you can see your data display on the screen, if right, click save current data if wrong, choose Re-measure.

5. Use App

In main screen, you can experience a lot of functions, such as Home management, trend checking, set target weight and so on, you can operate each function according to instruction in app.