**Escali®**

Measuring What's Next™

**Model: US200L**

**User's Guide**

**Easy Read  
Body Scale**

[www.escali.com](http://www.escali.com)

Thank you for choosing the Easy Read Body Scale by Escali. This precision state-of-the-art measuring instrument is most accurate when weighing on a flat hard surface. This scale will not function correctly on carpet. The scale's electronic sensors are sensitive; avoid dropping or jarring the scale and store it where it will be protected from impact.

## **Operating instructions**

1. Insert four AAA batteries (included) into the battery compartment under the scale.
2. Select your desired unit of measure by pressing the UNIT button on the top edge of the scale.
3. Place the scale on a solid and flat surface.
4. To activate simply stand directly onto the scale platform.
5. The scale will turn off automatically when finished.

## **Additional Information**

1. Make sure your scale is at room temperature, and on a stable horizontal surface free of vibration.
2. This scale is not waterproof. Do not expose to moisture.
3. As this scale is a precision instrument, it must be handled with extreme care. Avoid rough treatment, shaking and vibration.
4. Do not disassemble this scale or the warranty will be void.
5. The display value may be affected by electromagnetic disturbances such as operating near a radio. If this occurs, move the device that is causing the disturbance and turn the scale off and restart.

## **Product Support**

24/7 Online support: <https://escali.zendesk.com>

Email support: [feedback@escali.com](mailto:feedback@escali.com)

Phone support: 1-952-469-1965

## **Product Warranty**

Your Escali scale is guaranteed to be free of defects in materials and workmanship and protected by a Lifetime Limited Warranty. Please visit our website for more information. [www.escali.com](http://www.escali.com)

*Escali's liability is limited to two (2) times the cost of the product.*

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