

# Adjustable Adult Toilet

## Instruction Manual



✓ Anti Side-turning ✓ Thoughtful Design ✓ Strong Load-bearing ✓ Easy to Clean

## •Product Components



Main body 1 unit



Backrest 1pc



Armrest 2pcs



Legs 4pcs



Tissue Basket 1pc




Solid Bottom Inner Drum for Indoor Use



Squatting Toilet with Hollow Inner Bucket

Warm Tips: Inner drum type and PU seat cushion according to the purchase of the package configuration.

## •Installation Steps



- 1 Put the toilet body upside down, insert the legs into the proper position by aligning them with the holes, check and make sure that the legs are firmly installed and parallel.



- 2 Align the backrest with the holes as shown in the picture and press down firmly, check and make sure the backrest is firmly installed.



- 3 Align the armrests with the holes and press down firmly. Install the other armrest in the same way, check and make sure that the armrests on both sides are firmly installed.



- 4 Remove the lid and seat, put in the inner bucket and fix the paper towel basket.