







FIND YOUR SHAPE	FIND YOUR DIAMETER
 <p>KIDNEY</p> <p>Cabinet door is hinged and does not attach to lazy susan</p>	
 <p>PIE-CUT</p> <p>Cabinet door is inset and must attach to lazy susan*</p> <p><small>* If you are replacing an existing lazy susan, the existing susan will be attached to the door</small></p>	
 <p>FULL CIRCLE OR D-SHAPED</p>	
<p>① FIND YOUR SHAPE ② FIND YOUR DIAMETER*</p> <p>③ COMPARE MEASUREMENTS TO THE AVAILABLE PRODUCT SIZES/DIMENSIONS TO VERIFY FIT</p>	
<p><small>*Measure D1 & D2. Subtract 2" from the smaller of the two measurements to determine the largest diameter of Lazy Susan that can fit into the cabinet. (Example: If D1 is 32" and D2 is 30", max diameter is 28").</small></p> <p>NOTE: This product can be added to new or existing cabinetry or can replace an older model.</p>	
<p>WILL IT FIT?</p>	