

## **NEW!** Manipulative Cleaning Tub Kit (SAN-TUB1)



When students return, shared manipulatives like Lego®, tan gram shapes, counters, etc. will need a fast and easy way to be cleaned after use.

This simple but useful tub kit includes two of our durable Large Open Tubs and a new mesh-bottomed tub for a complete washing and drying set.

**Features:**

- 1 Mesh-bottomed tub
- 2 Large Open Tubs - 1 for washing and 1 for drying
- Handles on both sides
- Made of 25% recycled plastic

**NEW!** Manipulative Cleaning Tub  
(SAN-CC4068)

Don't need the whole kit? The Clear Mesh-bottomed Tub is available. Use on its own or with existing Large Open Tubs already in the classroom.


**Dimension of one tub**  
6" H x 12 1/2" W x 15 3/4" D

**Warranty**  
Lifetime




Kit shown (SAN-TUB1)

Our tubs have been tested with Caviwipes (medical grade disinfection wipes) and bleach diluted to 2:10




**How the kit works:**



**Step 1**  
Stack clear mesh tub inside the blue tub, add soap, water and manipulatives.



**Step 2**  
Once clean, drain, rinse with water in sink and drain.



**Step 3**  
Set on second blue tub to dry.