



Contemporary Copper Pull - 375

Antimicrobial

CuSalus® product - see below for advantages and benefits.

Product № 37596195AMC

Center to Center 96 mm

Length - Overall Dimensions 106 mm

SUGGESTED PRODUCTS



**Contemporary Copper Edge
Pull - 70 & 71**

Product № 7092195AMC
Center to Center: 2 in



**Contemporary Copper Pull -
8701**

Product № 870160195AMC



**Contemporary Copper Edge
Pull - 70 & 71**

Product № 7081195AMC
Center to Center: 1 in



**Contemporary Copper Pull -
75**

Product № 7534195AMC
Center to Center: 4 in



**Contemporary Copper Pull -
222**

Product № 22218195AMC
Screw/Nail: 1/4-20 (included)



**Contemporary Copper Pull -
222**

Product № 22228195AMC
Screw/Nail: Wood Screw
(Included)




**Contemporary Copper Pull -
2461**

Product № 246110195AMC



**Contemporary Copper Push
Plate - 94**

Product № 940816195AMC

 Contemporary Copper Pull - 75 Product № 7513195AMC Center to Center: 3 in	 Contemporary Copper Pull - 75 Product № 75235195AMC Center to Center: 3 1/2 in	 Contemporary Copper Pull - 75 Product № 7556195AMC Center to Center: 6 in	 Contemporary Copper Pull - 222 2 products
 Contemporary Copper Edge Pull - 70 & 71 2 products	 Breakaway Machine Screw, Large Truss Head, Combined Quadrex Slot Drive, 8-32, Type B Point Product № 2905 Packaging format: Bag of 50 units	 Breakaway Machine Screw, Large Truss Head, Combined Quadrex Slot Drive, 8-32, Type B Point Product № 2900 Packaging format: Bag of 100 units	

DISCLAIMER

The use of a CuVerro Copper Alloy surface is a supplement, not a substitute, for standard infection control practices. Users must continue to follow all current infection control practices, including those practices related to cleaning and disinfection of environmental surfaces. The Copper Alloy surface material has been shown to reduce microbial contamination, but it does not necessarily prevent cross-contamination.